



# NEVESKIN FACE

× ARTEMIS.

## Spec Sheet

### Key Features

- No downtime
- Easy to assemble
- Automatic feedback of client's temperature
- Minimal consumable costs
- Compact and portable
- Minimal noise
- Touch screen control
- 8 preset protocols

### Protocol Types:

- The Standard and Extended protocols: to brighten, contour, and lift.
- The Collagen Boost protocol: to boost and promote collagen and elastin production in the skin.
- The Regeneration protocol: to tighten and brighten.
- The Ultra Lift protocol: to lift, firm, and reduce signs of aging.
- The Calming protocol: to reduce inflammation and/or acne.
- The Age-Defying Upgrade protocol: to be used as an add-on to specifically target fine lines, wrinkles, and signs of aging.
- The Ultimate Eye Upgrade protocol: to be used as an add-on to specifically target the area around the eyes to brighten, de-puff, firm, and tighten the area around the eyes.

### Technical Specifications:

- Acoustic signal at the end of the session
- Voltage 100-240 VAC
- Frequency 50/60 Hz
- Absorption 200VA
- Output values: from 0 ° C to +41 ° C
- Protection class according to EN 60601-1: I B.
- Embedded software that does not require computer networks for its operation
- Software Class: B
- Dimensions- W10.86 inch, D14.57 inch, H12.40 inch,
- Weight- 4.309lbs